

HST-M2 Metallographic Specimen Pre-grinding Machine



Application:

HST-M2 type metallographic sample manual pre-grinding machine has smooth rotation, low noise, easy operation and high work efficiency. Suitable for metallographic laboratories of factories and scientific research institutions.

Feature:

- One machine can be used for multiple purposes, easily completing rough polishing and fine polishing of metallographic samples.
- Speed control system optimized for grinding and polishing machines
- Easy-to-operate switches and knobs for easy operation
- Finely ground and surface-treated work disc to ensure smooth sample surface
- The integrally formed high-strength fiberglass material shell is sturdy and durable and will never rust.

Main Technical Parameters:


Model	HST-M2
Polishing disc diameter	φ203mm
Rotating speed	450r/min
Net weight	40kg
Motor	370W
Dimension	712**660*330mm
Gross weight	58KG
Package size	758*776*421mm

Packing list:

Name	Specification	Qty
Metallographic sample polishing machine	HST-M2	1
Product Manual		1
Product certification		1
Packing List		1
Metallographic sand paper	Φ200	2
Grinding disc	Φ230	2
Water inlet pipe		1
Drain pipe		1
Water retaining ring		1

Metallographic consumables

Focus on metallographic equipment and consumables research and development of professional brand

			
Diamond cutting piece	Resin cutting sheet	Metallographic cutting Coolant	CBN cutting piece
			
CM-7 Cold mounting material	Cold/hot mosaic powder	Die for cold Mounting	Mounting material
			
Diamond Grinding Disk	Metallographic sandpaper	MEGULAR	Diamond paste
			
Polishing powder	Polishing lubricant coolant	Metallographic polishing cloth	Polishing liquid