

HST-XQ-2B Manual Metallographic Specimen Mounting Press



Application:

The XQ-2B mounting machine is used to mount metallographic and petrographic samples that are tiny, difficult to hold or irregular before grinding and polishing.

Features:

- One machine with multiple functions, easily completes pre-grinding, rough grinding, fine grinding, rough polishing and fine polishing of metallographic samples
- Stepless speed regulation, speed regulation control system optimized for grinding and polishing machines
- Easy-to-operate switches and knobs for easy operation
- Finely ground and surface-treated work disc to ensure smooth sample surface
- Integrated high-strength ABS material shell, sturdy and durable.

Technical Parameters:


Model	XQ-2B
The Sample Diameter	22mm OR 30mm OR 45mm (Please Choose One Sample Diameter)
Temperature control range	Room temperature-300℃
Timer Range	0-30min
Total power	≤600W
Dimension	265*260*425mm
Package Size	460*370*550mm
Net Weight	32kg
Gross Weight	40kg

Packing List:

Metallographic Specimen Mounting Press	XQ-2B	Set	1
Plastic spoon		Piece	1
Manual		Piece	1
Certification		Piece	1
Fuse		Piece	1
Thermocouple		Piece	1

Metallographic consumables

Focus on metallographic equipment and consumables research and development of professional brand

			
Diamond cutting piece	Resin cutting sheet	Metallographic cutting Coolant	CBN cutting piece
			
CM-7 Cold mounting material	Cold/hot mosaic powder	Die for cold Mounting	Mounting material
			
Diamond Grinding Disk	Metallographic sandpaper	MEGULAR	Diamond paste
			
Polishing powder	Polishing lubricant coolant	Metallographic polishing cloth	Polishing liquid