

HST-DL Cryogenic cooling circulation pump



Product Introduction

Widely used in universities, colleges, and enterprises in the fields of biology, chemical engineering, chemistry, pharmaceuticals, life sciences, and physical property testing. It provides users with a controlled, uniform, and constant temperature source for conducting constant temperature experiments or tests on experimental samples or manufactured products. It can also serve as a direct heating or cooling, or auxiliary heating or cooling heat source or cold source.

Technical Specifications:

Model	Temperature Range (℃)	Fluctuation (±℃)	Flow rateL/min	Maximum Head (m)	Volume	Cooling Capacity (W)
-------	-----------------------	------------------	----------------	------------------	--------	----------------------

DL-1005	-10-Room Temperature	±0.5	13	3	5	350
DL-2005	-20-Room Temperature	±0.5			5	180
DL-3005	-30-Room Temperature	±0.5			5	150
DL-4005	-40-Room Temperature	±0.5			5	200
DL-6005	-60-Room Temperature	±2			5	1910 ~ 300
DL-7005	-70-Room Temperature	±2			5	2020 ~ 250
DL-8005	-80-Room Temperature	±2			5	2100 ~ 210
DL-1010	-10-Room Temperature	±0.5			10	500
DL-1510	-15-Room Temperature	±0.5			10	600
DL-2010	-20-Room Temperature	±0.5			10	1150
DL-3010	-30-Room Temperature	±0.5			10	300
DL-4010	-40-Room Temperature	±0.5			10	400
DL-6010	-60-Room Temperature	±2			10	2750 ~ 300
DL-7010	-70-Room Temperature	±2	13	3	10	2870 ~ 280
DL-8010	-80-Room Temperature	±2			10	2900 ~ 230
DL-1015	-10-Room Temperature	±0.5			15	800
DL-2015	-20-Room Temperature	±0.5			15	700
DL-3015	-30-Room Temperature	±0.5			15	1000
DL-4015	-40-Room Temperature	±0.5			15	600
DL-6015	-60-Room Temperature	±2			15	2980 ~ 570
DL-7015	-70-Room Temperature	±2			15	3250 ~ 480
DL-8015	-80-Room Temperature	±2			15	3500 ~ 320
DL-1020	-10-Room Temperature	±0.5			20	1000
DL-2020	-20-Room Temperature	±0.5			20	1300
DL-3020	-30-Room Temperature	±0.5			20	1600
DL-4020	-40-Room Temperature	±0.5	13	3	20	800
DL-6020	-60-Room Temperature	±2			20	4120 ~ 510
DL-7020	-70-Room Temperature	±2			20	4380 ~ 420
DL-8020	-80-Room Temperature	±2			20	4530 ~ 360
DL-1030	-10-Room Temperature	±0.5			30	2300
DL-2030	-20-Room Temperature	±0.5			30	1800
DL-3030	-30-Room Temperature	±0.5			30	1000
DL-4030	-40-Room Temperature	±0.5			30	1000
DL-6030	-60-Room Temperature	±2			30	4330 ~ 680
DL-7030	-70-Room Temperature	±2			30	4520 ~ 570
DL-8030	-80-Room Temperature	±2			30	4679 ~ 440
DL-1050	-10-Room Temperature	±0.5			50	2500
DL-2050	-20-Room Temperature	±0.5			50	4000
DL-3050	-30-Room Temperature	±0.5			50	2800
DL-4050	-40-Room Temperature	±0.5			50	1500

DL-6050	-60-Room Temperature	±2			50	10570 ~ 1600
DL-7050	-70-Room Temperature	±2			50	10800 ~ 1480
DL-8050	-80-Room Temperature	±2			50	11200 ~ 1260

HST GROUP Co., LTD.

NO. 4915, West Jingshi Road, Jinan 250012, PR CHINA.

Tel: +86-15910081986

www.hssdgroup.com E-mail: admin@hssdtest.com