

## HST-GX High temperature circulator



### Product Introduction

Widely used in universities, colleges, and enterprises in the fields of biology, chemical engineering, chemistry, pharmaceuticals, life sciences, and physical property testing. It provides users with a controlled, uniform, and constant temperature source for conducting constant temperature experiments or tests on experimental samples or manufactured products. It can also serve as a direct heating or cooling, or auxiliary heating or cooling heat source or cold source.

### Technical Specifications:

| Model   | Temperature Range (°C) | Volume (L) | Heating power (kw) | Pump flow rate (L/min) | Fluctuation (±°C) | Resolution (°C) | Head (m) | External Dimensions (mm) |
|---------|------------------------|------------|--------------------|------------------------|-------------------|-----------------|----------|--------------------------|
| GX-2005 | RT-300                 | 5          | 1.3                | 13                     | 0.1               | 0.01            | 3        | 340*380*440              |
| GX-2010 | RT-300                 | 10         | 1.7                | 13                     | 0.1               | 0.01            |          | 340*380*490              |

|          |        |    |     |    |     |      |  |             |
|----------|--------|----|-----|----|-----|------|--|-------------|
| GX-2015  | RT-300 | 15 | 3.0 | 13 | 0.1 | 0.01 |  | 340*380*490 |
| GX-2020  | RT-300 | 15 | 3.5 | 13 | 0.1 | 0.01 |  | 340*380*490 |
| GX-2030  | RT-300 | 15 | 3.8 | 13 | 0.1 | 0.01 |  | 340*380*490 |
| GX-2050  | RT-300 | 20 | 4.4 | 13 | 0.1 | 0.01 |  | 340*380*620 |
| GX-20100 | RT-300 | 50 | 6.5 | 13 | 0.1 | 0.01 |  |             |

**HST GROUP Co., LTD.**

NO. 4915, West Jingshi Road, Jinan 250012, PR CHINA.

Tel: +86-15910081986

[www.hssdgroup.com](http://www.hssdgroup.com) E-mail: [admin@hssdtest.com](mailto:admin@hssdtest.com)