

Lever rubber tensile Grip HK-79



2KN,

Application:

Used for tensile testing of elastomers such as rubber and silica gel

Jaw 25*45mm ,Thickness 18mm

The tooth shape of the clip is mesh tooth (without hardening).

HST GROUP Co., LTD.

Jinan Hensgrand Instrument Co., Ltd.

NO. 4915, West Jingshi Road, Jinan 250012, PR CHINA.

Phone/WhatsApp: +86-159 1008 1986

www.hssdgroup.com - e-mail: admin@hssdtest.com