

HST-HPT3 Holding Power Tester



Product Uses

is manufactured in accordance with Method A of the national standard GB/T4851-2014 and is suitable for tack testing of various adhesive products, such as pressure-sensitive tapes, medical patches, self-adhesive labels, protective films , etc.

Technical parameters

1. Standard pressure roller: 2000g $\pm$ 50g
2. Weights: 1000 $\pm$ 1 g (including tripod weight) (200g optional)
3. Test plate: 125 (L) $\times$ 50 (B) $\times$ 1.5 (D) mm (60*40*1.5 mm optional)
4. Timing range : 0-99h 59min59s
5. Workstations: 3 pieces
6. External dimensions: 320 (L) $\times$ 210 (B) $\times$ 550 (H) mm
7. Control method: Microcomputer control; individual workstations can be switched on and off independently.

- 8. Display method: LCD monitor
- 9. Net weight: 12 kg
- 10 **】** Electricity Source: AC 220V , 50Hz

HST GROUP Co., LTD.

Jinan Hensgrand Instrument Co., Ltd.

NO. 4915, West Jingshi Road, Jinan 250012, PR CHINA.

Phone/WhatsApp: +86-159 1008 1986

www.hssdgroup.com - e-mail: admin@hssdtest.com